



OMI Mask pulsed electromagnetic therapy face mask



Portable PEMF Therapy Device for skin rejuvenation and increasing microcirculation in the face.

Description

Battery operated pulsed electromagnetic device specifically designed for the treatment of the face. The device runs on 4 AAA batteries. It contains 10 copper coils for the balanced distribution of the magnetic field on the face. The field intensity is 190 microT. Magnetic wave form is sinus.

Size: 10.24" x 12.20" (26 x 31 cm)

Recommended for

The OMI Mask is recommended literally for everybody over 40 wishing to preserve healthy shining look of the face and to avoid or reduce fine lines, wrinkles and to gain natural glow by increased micro-circulation.

Claims

The product provides electromagnetic massage on the targeted area. The pulsing magnetic field promotes the body's own regeneration processes and helps

- ✓ to prevent and treat dermatological problems,
- ✓ to increase microcirculation,
- ✓ to reduce inflammation,
- ✓ to reduce fine lines and wrinkles,
- ✓ to increase absorption of skin-care products.

Advantages

Easy to carry, for both professional and home use. Completely safe, not painful to use, and there are no known side-effects.

Instructions for use

The use of the product is very easy, just push the ON/OFF switch to turn on the device, and then place it onto the desired area you wish to treat.

How it works

Pulsed electromagnetic fields produced by the device induce electrical current inside the body that helps to restore the blocked ion flow and stimulates the body's own regeneration processes.

As a result of using the device you will attain healthy looking skin with improved blood flow.

Packaging contains

Pulsed electromagnetic device, operating unit, user manual.

Documentation and others

Declaration of Conformity (LVD, EMC)

Wholesale Inquires

All level of wholesale inquiries are welcome. To discuss further details, please [contact us](#).

Payment Terms: on delivery

Country of Origin: European Community

Pick-up Point: Budapest

Delivery Term, Currency and Price: [ask for quote](#)

Minimum Order Quantity: [ask for quote](#)

Package Details

Product Dimensions (cm): 29 x 26.5 x 4.5 cm

Product Net Weight (kg): 0.27 kg

Product Gross Weight (kg): 0.40 kg

Master Carton Quantity (pcs): 38 pcs

Master Carton Dimensions (cm): 64 x 60 x 46 cm

Master Carton Gross Weight (kg): 17.2 kg

The product is not a medical device; it is not intended to diagnose, treat, cure or prevent any disease.

product category: wellness product

website: www.oxfordmedicals.com/pemf-devices.html

